

GeneTrue™ (CP4 EPSPS) Test Kit

Instructions For Use

Format: Strip


Specimen: Extract of Soybean, Corn,
Cotton and other Crops

Catalog Number: A07-05-413



INTENDED USE

Artron GeneTrue™ Kit for Roundup Ready Crop is intended to test the presence of CP4 EPSPS protein at a level typically expressed in genetically modified (GM) Roundup Ready Crop (Soybeans, Corn, Cotton, etc) to assist the determination of Roundup Ready Crops. The sensitivity of the test is 1 ng/ml, or 0.1% (i.e. one Roundup Ready seed in 1000 conventional kernels).



This assay provides only a preliminary result. Professional judgment should be sought to further evaluate the result of the test, particularly in evaluating a preliminary positive result.


SUMMARY AND PRINCIPLE OF THE ASSAY

CP4 EPSPS is the protein introduced by genetic manipulation that is expressed in glyphosate-tolerant soybeans, which are being developed to provide new weed-control options for farmers. Expression of this protein in plants imparts high levels of glyphosate tolerance. The safety of CP4 EPSPS was ascertained by evaluating both physical and functional characteristics. CP4 EPSPS degrades readily in simulated gastric and intestinal fluids, suggesting that this protein will be degraded in the mammalian digestive tract upon ingestion as a component of food or feed.

Artron GeneTrue™ Kit for Roundup Ready Crop is an antigen-capture immunochromatographic assay, which detects the presence of CP4 EPSPS in extract of soybeans or other crop sample. Monoclonal antibodies specifically against CP4 EPSPS are 1) conjugated with colloidal gold and deposited on the conjugate pad, and 2) immobilized on the test line of the nitrocellulose membrane. Artron GeneTrue™ Strip has an absorbent pad at each end. The protective tape with the arrow indicates the end of the strip to insert into the samples.

In order to detect the CP4 EPSPS protein expressed by Roundup Ready Crops, the sample must first be extracted to solubilize the protein. After the strip inserted into the sample, the gold-antibody conjugate is rehydrated and the CP4 EPSPS, if any in samples, interacts with the gold conjugated antibodies. The antigen-antibody-gold complex will migrate toward the test window until the Test Zone (T) where it will be captured by immobilized antibodies, forming a visible pink line (Test band), indicating a positive result. If CP4 EPSPS is absent in the sample, no pink line will appear in the Test Zone (T), indicating a negative result.

To serve as an internal process control, a control line should always appear at Control Zone (C) after the test is completed. Absence of a pink control line in the Control Zone is an indication of an invalid result.



PACKAGE CONTENTS

- 100 Artron GeneTrue™ Strips packed in two canisters (50 strips/canister, 2 canisters/box)
- 1 Instructions for use
- 100 transfer pipettes and 100 reaction vials

MATERIALS REQUIRED (BUT NOT PROVIDED)

- Clean, dry collection container (plastic or glass).
- Clock or timer.
- Sample grinder.
- Pestle (for leaf samples only)
- Tap water

WARNINGS AND PRECAUTIONS

- Do not reuse.
- Do not use if the canister or its packaging is compromised.
- Do not use after the expiration date shown on the packaging.
- Do not mix and interchange different specimens.
- Wash hands thoroughly after finishing the tests.
- Do not eat, drink or smoke in the area where the specimens or kits are being handled.
- Keep out of children's reach.

SPECIMEN PREPARATION

1. Determine Number and Size of Sub-samples

- For sampling plan, please read literatures in the Reference Sections. Collect a composite sample according to USDA/ GIPSA instructions found in the reference literature.
- To select an appropriate sample size, determine the purity standard and the degree of confidence required. Confidence level means the statistical probability that the true Roundup Ready level in the lot is below the selected purity standard. Table 1 provides a guideline for determining the number of beans in each sub-sample that are necessary to provide effective screening for different GM concentrations at the 95% and 99% confidence levels.

Table 1 Soybeans Number of beans per sub-sample at different GM concentrations

1-1: Testing 2 sub-samples one or both sub-samples being negative (-/-) or (-/+)

Confidence Level (%)	Roundup Ready Screening Level Soybeans						
	5%	4%	3%	2%	1%	0.5%	0.4%
95%	72	90	120	182	366	734	915
99%	104	130	175	262	527	1060	N/A

1-2 Testing 3 sub-samples two or all three sub-samples being negative (-/-) or (-/-+)

Confidence Level (%)	Roundup Ready Screening Level Soybeans						
	5%	4%	3%	2%	1%	0.5%	0.4%
95%	39	49	66	100	200	400	1000
99%	55	70	93	141	281	565	N/A

Note: Screening at the 0.1% Roundup Ready concentration level, with 95% confidence, would require testing 3 sub-samples of 1000 beans with all 3 sub-samples negative (-/-).

2. Determine Sub-Sample Weight, and Grind Times

- Determine average weight of individual grain to be tested (weigh 100 beans divide by 100).
- Calculate the weight of the number of grains to be tested (Number of grains x Average Weight/Grain). Table 2 lists the guidelines for grinding time according to sample weight.

Table 2: Grinding time and sample weight.

Commodity	Sample Weight (g)	Grind Time (Sec)
Soybeans	9-15	20
	16-40	20
	41-60	30
	61-150	45

3. Determine the Volume of Water to be Used in Extraction

The test is suitable for testing both leaf tissue and seed of Soybean, corn, cotton and other crops.

Leaves, or seeds must be grinded and diluted with water (Tap water could be used). For best result, different samples should be diluted with water in different ratio (see table below). Let the paste set for at least 45 seconds before testing with the test kits.

To extract seed:

A single seed sample can be grinded and then extracted with water in a tube. The sample preparation is a crucial procedure for the proper function of the test. The ratio of water in milliliters (ml) to the weight of the seed sample (g) ranges from 1:5 to 1:25 depend on the different crop.

To extract leaf tissue:

Take a 1-inch x 1-inch leaf punches from the sample leaf tissue; push the leaf punches into the tube by a pestle. Then squeeze 10 drops of water into tube and grind the tissue using the pestle.

Table 3: Sample Weight and Water Ratio

	Weight of Sample (g)	Volume of Water (ml)	Ratio weight (g):water (ml)
Soybean seed	20	100	1:5
Corn seed	20	100	1:5
Cotton seed	5	100	1:20
Single Corn Leaf	5	100	1:20
Multiple (composite) Leaves	10	100	1:10

4. Prepare the Sample

- Weigh sample into the appropriate size glass or plastic jar.
- The sample can be grinded either dry or in water depend on the grinder available.
- Grind sample with a blender (on high speed) for specified grinding time or until all whole grains are broken into powder.
- Add the volume of tap water calculated according to Table 3 above.
- Cap and shake jar vigorously until the entire sample is wet (20-30 seconds, depending on the number of grains). Sample will begin to settle immediately and liquid can be drawn off at that time.
- Transfer of about 1 ml of the liquid of the sample to the reaction vial with the transfer pipette. (Draw up enough liquid portion from above the settled sample to fill the long narrow tip of the transfer pipette and up to the line at the top of the flared portion of the pipette bulb). Avoid pulling up particles.




Important: To prevent cross-contamination thoroughly clean blender parts and jars of dust and residue prior to preparation of a second sample. Use a new transfer pipette and reaction vial for each sample.

TEST PROCEDURES

1

Allow refrigerated canisters to come to room temperature before opening. Remove the canister from the foil pouch by tearing at the notch. Remove the Artron GeneTrue™ Strips to be used. Reseal the canister immediately.



2

Hold the strip at the colored end, and immerse the strip into the specimen in the reaction vial with the arrow end pointing towards the specimen. Do not immerse past the MAX line.



3

You could leave the strip in the reaction vial until the test is completed, or take the strip out after a minimum of 10 sec. Lay the strip (MAX side facing up) flat on a clean, dry, non-absorbent surface.



4


Read the result in 5 minutes. Read results as shown under interpretation of Results.

NOTE: Specimens with high concentrations of CP4 EPSPS may produce positive result in as little as 1 minute. Confirm negatives within 10 minutes



**DO NOT INTERPRET RESULTS
AFTER 30 MINUTES**

RESULT INTERPRETATIONS



Negative

A pink band appears only at the control region, indicating absence of CP4 EPSPS or the concentration of the CP4EPSPS in the sample is below detection level (1 ng/ml).

Positive

Two pink bands appear at the control and test regions, indicating a positive result for CP4 EPSPS.

Invalid

No visible band appears at the control region. Repeat with a new test kit. If the test still fails, please contact the distributor with the lot number.

QUALITY CONTROL

Although the testing device contains an internal quality control (pink colored line in the control region), good laboratory practice recommends the daily use of an outside control to ensure proper testing device performance. Quality control samples should be tested according to the standard quality control requirements established by your laboratory.

STORAGE AND STABILITY

- Test device in the canister should be stored at 2-30°C. Do not freeze the test device.
- The test device should be kept away from direct sunlight, moisture and heat.

LIMITATIONS

- Humidity and temperature can adversely affect results.
- The instructions for the use of the test should be followed during testing procedures.
- There is always a possibility that false results will occur due to the presence of interfering substances in the specimen or factors beyond the control of the manufacturer, such as technical or procedural errors associated with the testing.
- Although the test demonstrates superior accuracy in detecting CP4 EPSPS, a low incidence of false results can occur. Therefore, other available tests are required in case of questionable results.

REFERENCES

1. <http://www.archive.gipsa.usda.gov/reference-library/handbooks/graininsp/grbook1/bk1.pdf>- USDA Grain Inspection Handbook, Book 1, Grain Sampling.
2. <http://www.archive.gipsa.usda.gov/biotech/sample2.htm>- Guidance document entitled Sampling for the Detection of Biotech Grains.
3. <http://www.archive.gipsa.usda.gov/biotech/sample1.htm>- Practical Application of Sampling for the Detection of Biotech Grains.
4. <http://www.archive.gipsa.usda.gov/biotech/samplingplan1.xls>- Sample Planner Excel Spreadsheet. The planner allows you to enter different assumptions in terms of sample size, number of samples, acceptable quality level and to determine the probability of accepting lots with given concentration levels.
5. Remund, K.M., Dixon, D.A., Wright D.L., Holden, L.R. "Statistical considerations in seed purity testing for transgenic traits", Seed Science Research, June 2001, Vol. 11 No.2, pp. 101-119.

FOR TECHNICAL ASSISTANCE OR QUESTIONS

Please contact us at:

Phone: 604-415-9757

Fax: 604-415-9795

Email: info@artronbio.com



LIMITED WARRANTY

Artron Laboratories Inc. ("Artron") warrants the products sold hereunder ("the Products") against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product's printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies Artron in writing of such defects during the warranty period, including an offer by the customer to return the Products to Artron for evaluation, Artron will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

NEITHER Artron NOR ITS PARTNERS MAKES ANY OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of Artron products appearing in Artron published catalogues and product literature are Artron' sole representations concerning the Products and warranty. No other statements or representations, written or oral, by Artron' employees, agents or representatives, except written statements signed by a duly authorized officer of Artron Laboratories Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

Artron does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by Artron, or against damages resulting from such non-Artron made products or components. Artron passes on to customer the warranty it received (if any) from the maker thereof of such non-Artron made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Artron.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of Artron shall be to repair or replace the defective Products in the manner and for the period provided above. Artron shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall Artron be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Artron with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

MANUFACTURER CONTACT INFORMATION



Artron Laboratories Inc.
3938 North Fraser Way
Burnaby, BC
V5J 5H6 Canada

Ph: 604.415.9757
Fax: 604.415.9795
www.artronlab.com
info@artronlab.com

Doc No. A07-05-413-CJ155EN

Version No: 15-01

Revised at August, 2015